




EMG 5000Q EMG/NCV/EP

**Accurate
Reliable
Easy To Use**

NR Sign Manufacturing 5000Q/NCV/EP is a full featured diagnostic system with a low noise high performance amplifier and state of the art software features. This system can be configured for 2 or 4 channels.

F Wave
H reflex
Blink Reflex
RNS (Repetitive Nerve Stimulation)
Motor NCV
Sensory NCV
Needle EMG
Visual Evoked Potential
Auditory Evoked Potential
Somatosensory Evoked Potential
CE, ISO 13485, CMDCAS, FDA QRS certified



**Accurate
Reliable
Easy to use**

**Model:5000Q
EMG/NCV/EP**




Protocol wizard:
This is a powerful tool for defining new tests, designing and customizing your system in your own way.

Name of nerves, muscles, number of signal traces, method of calculations, shape of result Table and final report and many other options can be defined by Protocol wizard.



- Microsoft® Word report generator can help you to generate your report in less than a minute.
- Graphical directory for muscle and nerve electrode placement can give you useful recommendations.
- There are many of applicable functions like:
 - Left and Right comparison
 - Super imposer
 - Pick detector
 - Needle Video and Audio reviewer
 - Signal annotation to help you to get more precise results.



- Microsoft® SQL Server Guarantee the integrity of data management this data can be accessed trough network.